



PARASOL MEDICAL, LLC

1602 Barclay Blvd • Buffalo Grove, IL 60089
Phone: 847-520-3000 • Fax: 847-520-4268



MICROBECARE™
Simplify for better care - MicrobeCare™

| MicrobeCare Cleaning Materials | | | | |
|---------------------------------------|---|-------------------|------------------------------------|--|
| Test Solution | Description | Tested Y/N | **Recommended Exposure Time | Active Ingredient |
| <i>Cidex</i> | <i>High Level Disinfectant Solution</i> | <i>N</i> | <i>12 minutes</i> | <i>Ortho-Phthalaldehyde (OPA) .55%</i> |
| Sani-Cloth | Germicidal Disposable Cloth | Y | 3 minutes | Dimethyl ethylbenzyl ammonium chloride .14%, dimethyl benzyl ammonium chloride .14% |
| A-456 II | Disinfectant Cleaner | Y | 10 minutes | Octyl decyl dimethyl ammonium chloride 6.510%, Dioctyl dimethyl ammonium chloride 2.604%, Didecyl dimethyl ammonium chloride 3.906%, Alkyldimethyl benzyl ammonium chloride 8.680% |
| Re-Juv-Nal | Disinfectant Detergent | Y | 10 minutes | Octyl decyl dimethyl ammonium chloride 1.62%, Dioctyl dimethyl ammonium chloride .814%, Didecyl dimethyl ammonium chloride .814%, Alkyldimethyl benzyl ammonium chloride 2.17% |
| Precise | Hospital Cleaner Disinfectant | Y | 10 minutes | 0-Phenglphenol .37% |
| Bleach | Household Cleaner | Y | N.A | Sodium Hypochlorite 100% |
| Super HDQ Neutral | One Step Germicidal Detergent | Y | 10 minutes | Didecyl dimethyl ammonium chloride 9.22%, n-Alkyl (C ₁₄ 50%, C ₁₂ 40%, C ₁₆ 10%) dimethyl benzyl ammonium chloride 6.14% |

** This list is a partial list of EPA approved disinfectant solutions commonly used in healthcare environments and approved to be used on MicrobeCare treated surfaces. If your institution uses a disinfectant not contained herein please contact us to inquire about its compatibility.*